

FREEZE-600ML FREEZING COLORLESS 600ML AG TERMOPASTY



New AG TERMOPASTY FREEZE-600ML Tested and cleaned. 24 months guarantee Dedicated courier Delivery across Europe even in 24h

TECHNICAL DATA

MANUFACTURER	AG TERMOPASTY
MODEL	FREEZE-600ML
CATEGORY	PLC subassemblies
WEIGHT [KG]	1.0
HEIGHT [CM]	10.0
LENGHT [CM]	33.0
WIDTH [CM]	10.0





TEST STAND



We have test stands developed by us, thanks to which we check whether the device is in working order. It is with these posts that we test the repaired device, for which we give 24 months warranty, and you can be sure that your FREEZE-600ML equipment has been repaired at a service center that provides the highest quality services

PRODUCTS FROM THE RGB ELEKTRONIKA OFFER:

HIGHEST QUALITY



assessment, based on which we select the best quality equipment to guarantee their trouble-free operation. Thanks to this, you can be sure that the AG TERMOPASTY FREEZE-600ML device has been prepared with the greatest care.

EXPRESS SHIPPING



Each product undergoes a specialized We send products from our warehouse the same day the order was placed. We also send a dedicated courier on request. You can be sure that your AG TERMOPASTY FREEZE-600ML will reach you as soon as possible.

24 MONTHS GUARANTEE



All used products are covered by our 24 months warranty! This means that the AG TERMOPASTY FREEZE-600ML device you purchased will be repaired by us under warranty for 2 years, and you can be calm about the smooth flow of your business.







COOPERATION WITH RGB ELECTRONICS IS:

QUICK SUPPORT



At RGB Elektronika, we understand that every minute counts in an emergency. Choose AG TERMOPASTY FREEZE-600ML from our offer and we will deliver them to you even the same day. Thanks to the express cooperation with our company, your company will be able to operate without downtime!

OWN SERVICE



We repair devices on our website in Wroclaw. We have the most modern equipment supported by experienced technicians. These independent activities help us to work faster and maintaining the highest quality of services. You can be calm about the efficiency of AG TERMOPASTY FREEZE-600ML.

REPAIR HISTORY



On request, we can attach a repair history to AG TERMOPASTY FREEZE-600ML. Thanks to this, every customer can see which elements have been replaced and which tests we have subjected to. At RGB Elektronika, we value honesty. We want our clients to trust us and be satisfied with the equipment purchased from us.

TRANSPORT SAFETY



You don't have to worry about damaging AG TERMOPASTY FREEZE-600ML during transport. We know how important it is to safely deliver a product to the customer. Thanks to many years of cooperation with UPS, we have developed innovative ways of securing devices in transport. We use three types of packaging for shipped devices: wooden box, instapak, dedicated pallet.

DO YOU KNOW WHAT TO DO WHEN A DAMAGED SHIPMENT WILL BE DELIVERED TO YOU?



Report damaged parcel to us within 7 days. You can do it by sending us an email message: biuro@rgbelektronika.pl or by calling our hotline: +48 71 325 15 05.



Take a picture of the damaged package and attach it to the email. Do not dispose of the packaging, invoice and shipping documents!



Don't worry about anything, we'll take care of everything. We will arrange a courier who will come to write you a damage report.



You can count on us - we will take care of everything for you!

We will inform you about each stage of the complaint procedure.

LET'S STAY IN CONTACT

RGB Elektronika spółka z ograniczoną odpowiedzialnością sp. k. Jana Długosza 2-6 51-162 Wrocław Poland



24/7 +48 71 750 09 78





